



STANWADE METAL PRODUCTS

the best made is Stanwade

To whom it may concern:

This letter applies to all makes and models of steel tanks that we have manufactured in any time period, including all:

- Single-wall AST & UST
- Double-wall AST & UST
- STI-P3
- ACT-100U
- Non-UL storage vessels

All steel tanks are compatible and suitable for use with all blends of fuels meeting ASTM standards, including ethanol blends from E10 to E100. All tanks are also compatible and suitable for use with all blends of biodiesel, from B2 to B100. Testing has been done proving compatibility of steel by several sources including Oak Ridge National Lab sponsored by DOE in collaboration with UL and NREL. SwRI, DNV and STI (through Battelle). To access test reports and other information on biofuels, see:

[STI .com website for test data and information on biofuel testing and steel compatibility](#)

As always, if switching product service, the tank system should be properly cleaned. See STI Recommended Practice, R111, "Storage Tank Maintenance" for further information.

If you have any questions about our products or about this statement please contact:

Timothy K. Woofter
President,
Stanwade Metal Products, Inc.
(800) 826-5243
6868 State Route 305
P.O. Box 10
Hartford, Ohio 44424



Timothy K. Woofter, President

6868 St. Rt. 305 • P.O. Box 10 • Hartford, OH 44424
330.772.2421 • 800.826.5243 • Fx 330.772.3307

www.tankstore.com





STANWADE METAL PRODUCTS

the best made is Stanwade

To whom it may concern:

This letter applies to all makes and models of steel tanks that we have manufactured in any time period, including all:

- Single-wall AST & UST
- Double-wall AST & UST
- STI-P3
- ACT-100U
- Non-UL storage vessels

All steel tanks are compatible and suitable for use with all blends of fuels meeting ASTM standards, including ethanol blends from E10 to E100. All tanks are also compatible and suitable for use with all blends of biodiesel, from B2 to B100. Testing has been done proving compatibility of steel by several sources including Oak Ridge National Lab sponsored by DOE in collaboration with UL and NREL. SwRI, DNV and STI (through Battelle). To access test reports and other information on biofuels, follow this link:

Biofuel data, information and links

As always, if switching product service, the tank system should be properly cleaned. See STI Recommended Practice, R111, "Storage Tank Maintenance" for further information.

If you have any questions about our products or about this statement please contact:

Timothy K. Woofter
President,
Stanwade Metal Products, Inc.
(800) 826-5243
6868 State Route 305
P.O. Box 10
Hartford, Ohio 44424



Timothy K. Woofter, President